

PHYSIOLASER 500

LASER THERAPY

THE PRODUCT IS INDICATED TO ACCELERATE THE RECOVERY OF INFLAMED TISSUES, REGENERATE CELLS AND IMPROVE MICROCIRCULATION. PHYSIOLASER 500 PROMOTES THE NATURAL HEALING PROCESSES.

TECHNICAL SPECIFICATIONS

Display

2,6", backlit

Laser source

Class 4 GaAlAs diode

Maximum power

0,5W ±20%

Compatible with a wavelength of

808nm

Issue

Pulsed and continuous

Target area dimensions

11 mm optic

Frequency

Up to 10,000Hz

Power supply


Mains and rechargeable batteries

Size and weight

100x160x35 mm - 433 g

STANDARD EQUIPMENT

- 1 Suitcase
- 1 PHYSIOLASER 500 device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



CE
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

GLOBUS
ITALIAN EXCELLENCE

PHYSIOLASER 500

LASER THERAPY

PROGRAM LIST

HEALTH - REHAB

Arthrosis

Ankle sprain

Knee pain

Edema Ankle sprain

Shoulder articular pain

Epicondylitis

Chronic pain Back pain

Neck pain

Tendinopathies

AESTHETIC - FITNESS

FACE AESTHETIC

Wrinkles

Acne

Skin rejuvenation

Skin - face lift

Expression wrinkles

Rosacea

BODY AESTHETIC

Breast stretch marks

Thighs, hips and gluteus edematous cellulite

Non-localized body fat and cellulite

Stretch marks

Scars